




OPHTHALMOLOGY INSTRUMENTS

www.praise.com.pk



Straight Tying Forceps




Straight shafts with 4.5mm tying platforms, 85mm long

Tips

- E-001** Standard
- E-002** Ultrafine

Mcpherson Tying Forceps



Angled shafts with 4.5mm tying platforms, 85mm long

Tips

Shafts

- | | | |
|--------------|-----------|------|
| E-003 | Standard | 12mm |
| E-004 | Ultrafine | 07mm |
| E-005 | Ultrafine | 10mm |
| E-006 | Ultrafine | 12mm |

Castroviejo Tying Forceps




5.0mm long smooth jaws for 8/0 to 11/0 suture. Shafts available straight or angled 45°, Castroviejo handle. Overall length 108mm

Shafts

- E-007** Straight
- E-008** Angled

Straight Tying Forceps




Straight shafts with 4.5mm tying platforms, 115mm long

Tips

- E-009** Standard
- E-010** Ultrafine

Mcpherson Tying Forceps



Angled 45° shafts with 4.5mm tying platforms 115mm long

Tips

- E-011** Standard
- E-012** Ultrafine

**McPherson Forceps
E-088**



**McPherson Forceps
E-089**



**Suture Tying Forceps
E-090**



**Suture Tying Forceps
E-091**



**Utrata Capsulorhexis Forceps
E-092**



**Utrata Capsulorhexis Forceps
E-093**



**Corneal Forceps 1x2 (0.12)
E-094**



**Corneal Forceps 1x2 (0.12)
E-095**



DeBekay Forceps
E-096
 5.5"




DeBekay Forceps
E-097
 5.5"



Serrated Curved Forceps
 with lock
E-098
 5.5"



Adson Tissue Forceps Titanium
 1.0mm tips / 4 3/4" (12cm)
E-099 TC dust
E-100 Tying Platform



Adson-Harris Forceps Titanium
 Fenestrated / ring-tip /
 TC dust / 4 3/4" (12cm)
E-101



Adson Tissue Forceps Titanium
 1 x 2 Teeth / 5 7/8" (15cm)
E-102



Brown Adson Tissue Forceps Titanium
 1.8 mm tips / 7x7 teeth / 4 3/4" (12cm)
E-103



Troutman Micro Hook Forceps
 Titanium with Tying Platform / TC dust / round Handle
E-104





Akahoshi Prechopper



Conventional smooth sharp pointed tips, 12mm angled curved shafts cross-action, maximum tip opening 3.0mm

E-119 120mm long



Akahoshi Phaco Prechopper



4.5mm long delicate jaws with conventional smooth sharp tips for cracking dividing and rotating the nucleus without sculpting, 120mm long

E-120 Shafts Straight
E-121 Angled

E-122 Nucleus Bisector/Vectis Forceps



Nucleus Trisector/Vectis Forceps



E-123 110mm long

CONJUNCTIVA FORCEPS



Fechtner Conjunctiva Forceps



Very delicate ring tips, with tying platforms

E-124 112mm long

CHALAZION FORCEPS

E-125 Chalazion Forceps



Oval 12 X 14mm, 98mm long



E-126
Round 8mm,
98mm long



E-127
Round 15mm,
98mm long




E-128
Oval 17 X 9mm,
98mm long



E-172 Castroviejo Needle Holder



Castroviejo Needle Holder



Fine 9mm straight jaws, wide serrated handle, sharp tips, 125mm long

	Lock
E-173	Without lock
E-174	With lock

Castroviejo Needle Holder



Fine 9mm straight jaws, wide serrated handle, sharp tips, 125mm long

	Lock
E-175	Without lock
E-176	With lock

SCISSORS

E-177 Mcpherson-Westcott Conjunctiva Scissors



Curved, blunt tips. Tip to Pivot length 10.5mm 105mm long


E-178 Conjunctiva Scissors



Gently curved blades. Tip to pivot length 12mm. Blunt tips, 115mm long



E-179 Osher Universal Scissors



Thin curved blades. Blunt tips, 21mm from pivot to tip, 115mm long

CORNEAL

Castroviejo Corneal Section Scissors

E-180



Right

E-181



Left

Blunt tips and lower blade 0.5mm longer than the upper blade, 13mm from pivot to tip, 105mm long

E-182 Castroviejo Keratoplasty Scissors



Universal with blunt angled blades. Blunt tips and lower blade longer than the upper blade, 16mm from pivot to tip, 105mm long

Castroviejo Corneal Section Scissors

E-183



Right

E-184



Left


Blunt tips and lower blade 0.5mm longer than the upper blade, 9mm from pivot to tip, 105mm long

E-185 Castroviejo Curved Corneal Scissors



Universal medium model with blunt tips. 16mm from pivot to tip, 116mm long

E-186 Castroviejo Corneal Scissors



Curved blunt tips, small blades, 105mm long

E-187 Castroviejo Corneal Scissors



16mm from pivot to tip, 116mm long

E-188 Castroviejo Corneal Scissors



Blunt tips, curved blades, 11.0mm from pivot to tip, 112mm long

E-189 Castroviejo Curved Scissors



Blunt tips, tip to pivot length 11.0mm, 105mm long

E-190 Castroviejo Straight Scissors



Blunt tips, tip to pivot length 11.0mm, 105mm long



Phaco Chopper



Phaco side cutter chooper, rounded and highly polished distal tip

Blade
E-219 Right hand
E-220 Left hand



Phaco Chopper



Phaco side cutter chooper, rounded and highly polished distal tip

Blade
E-221 Right hand
E-222 Left hand


E-223

Minami "m" Hook



Chopper, Splitter, Dissector, Spatula, Rotator and Manipulator in one

Gillum Nucleus Splitter



E-224

E-225

E-226

Claw Phaco Chopper



Phaco Chopper



Ambidextrous nucleus chopper, tip 1.8mm long
 Rounded and highly polished distal tip

Blade
E-227 Straight
E-228 Cuved



Phaco Chopper



Phaco side cutter chopper, rounded and highly polished distal tip
 Tip 1.8mm long

Blade
E-229 Right hand
E-230 Left hand

E-231 Haefliger Phaco Cleaver



Chopper and manipulator in one

E-232 Phaco Chopper



Special inner blade is offset 45°, rounded distal tip is highly polished, unique elongated distal segment

Phaco Chopper



Rounded and highly polished distal tip. Tip 1.8mm long

Blade
E-233 Right hand
E-234 Left hand

E-235 Phaco Chopper





Paddle manipulator / 60° Phaco chopper



E-236 Phaco Chopper



E-237 Schocket Double Ended Scleral Depressor



E-238 Phaco Chopper



145mm long






E-291 Akahoshi Nucleus Stainer/manipulator

TITANIUM BLUE COATED



E-292 Intracapsular Manipulator



E-293 Intracapsular Manipulator




RETINA HOOKS AND SPATULAE

E-294 Retina Hook




Triangular 0.5mm tip, angled 130°

E-295 Membrane Spatula



Rounded tip spatula

E-296 Retina Spatula



Tip angled 110°

E-297 Retina Hook



E-298 Kansas Lens Loop

TITANIUM BLUE COATED



3mm wide. Serrated. For small incision nucleus fragment removal

E-299 Kansas Nucleus Bisector



For manual small incision surgery

E-300 Kansas Nucleus Trisector



For manual small incision surgery

E-301 Kansas Nucleus Vectis



Phaco Chopper



Shafts
E-302 Straight
E-303 Angled

E-304 Nealiann Nucleus Divider



Ambidextrous nucleus chopper

Nichamin Phaco Chopper



Ambidextrous nucleus chopper, tip 1.1mm long.
Rounded and highly polished distal tip


Blade
E-305 Straight
E-306 Cuvéd

E-307 Spatulated LASIK Cannula



This 26G curved cannula has an end opening for irrigating the interface. Flattened spatulate tip to easily locate and elevate the flap edge.

E-308 Gimbel LASIK Fountain Cannula



Smooth bullet shaped tip single centerport for directing fluid toward the superior hinge, 25G

E-309 Buratto LASIK Cannula

25G



One front and two side ports. Designed for equal flow in three directions

E-311 Banaji LASIK Cannula



Bullet shaped tip. 6 ports for multi-directional irrigation, 25G

E-312 Vidaurri LASIK Cannula



Double-armed 8 irrigating ports 25 G

CANNULA-GENERAL


Posterior Capsule Polisher

E-313 25 gauge
E-314 23 gauge



Titanium impregnated on posterior surface of olive tip. Front opening. Gently curved

E-315 Air Injection Cannula




Angled 45° shaft, 5mm from tip to bend with polished finish, 27G

E-316 Nichamin Hydrodissection Cannula




Curved cannula. 25G

E-317 Lacrimal Cannula



23G. Cannula with blunt tip. Straight shaft.

E-318 Lacrimal Cannula




23G. Cannula with reinforced shaft. Straight 10mm tip.

E-319 Lacrimal Cannula



23G. Cannula with reinforced shaft. Curved 10mm tip.

E-320 West Lacrimal Cannula



23G. Cannula with blunt tip and bottom port. Gently curved shaft



Excellence in Quality Instruments

If you need full catalogue please contact us

Praise Corporation

Address : Muslim Town, Airport Road Kotli Behram Sialkot 51310
Pakistan.

Tel : +92-52-4273725

+92-52-4273723

+92-52-4264934

Cell / Whatsapp : +92 301 8606300

Fax : +92-52-4262310

Email : info@praise.com.pk

Muhammad Afzal Mughal (Director)

Email : afzal@praise.com.pk

